# CD63/CD63 Exosome ELISA Kit, Human

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---

#### **【 I – 1 】 Background and Measurement Principal**

Exosomes are membrane vesicles, which are secreted from many cells and have a diameter of about 30 nm ~ 200 nm. They are present in many bodily fluids, including saliva and blood <sup>1) ~ 5)</sup>. It is also known that exosomes are secreted into the mammalian cell culture medium in vitro. Exosomes are enveloped with a lipid bilayer membrane, so like cells, membrane proteins exist at their surface, and they encase various molecular constituents, including proteins and microRNAs. When the target cells receive the exosomes, the proteins or microRNAs inside will fulfill their function, and in this way exosomes are thought to play a role of cell-to-cell communication <sup>6)</sup>. One of the exosomes' structural characteristics is the tetraspanin family which is located on their surface. CD9 and CD63 are member molecules of this family, and known to be exosome surface markers <sup>7)</sup>.

The conventional way to quantify exosomes is indirect. For example, it measures encapsulated protein abundance or performs Nanoparticle Tracking Analysis to determine the size distribution profile of small particles <sup>8)</sup>. The drawback for these methods is that they require exosome purification utilizing ultracentrifugation, for instance. The direct methods to quantitate exosomes in body fluids or culture supernatant is extremely limited, and there were no common methods available until now. The CD63/CD63 Exosome ELISA Kit is a Sandwich ELISA kit, which utilizes high-performance anti-CD63 antibodies. This product detects exosome markers, CD63 molecules that are located on the exosome surface in the body fluids or cell culture supernatant.

## 【 I – 2 】 Features

- Directly quantitate exosomes in human blood samples or cell culture supernatant.
- No special equipment is required. Standard microplate reader capable of reading at 450nm will do the job.
- Utilize CD63 immobilized beads of 200 nm diameter(CD63 Standard Beads), instead of unstable/hard to store exosome itself, to implement stability and reproducibility.
- Normalization with CD63 standard beads similar to exosomal structures (CD63 Standard Beads) enable to relative quantitate each samples.
- Capture exosomes with solid phase anti-CD63 antibody, then detect using HRP conjugated anti-CD63 antibody.

## 【 I – 3 】 Kit Principle

An anti-CD63 antibody was immobilized and placed on the plate. First, samples were added onto the plate to capture exosomes by the anti-CD63 antibody. Next, HRP conjugated anti-CD63 antibody will be added to react exosome surface antigen, CD63. Finally, substrate will be added, then measure the coloring by the plate reader to quantitate sample exosomes.

## 【 I – 4 】 Kit Component

Storage temperature : 2 ~ 8 °C

	Reagent	Volume	Quantity
1	Anti-CD63 Antibody Immobilized Plat	96well (8well x 12 strips)	1 plate
2	CD63 Standard Beads	200μL	1tube*
3	Assay Buffer	25mL	1vial
4	Washing Buffer (10X)	25mL	1vial
5	HRP Conjugated Anti-CD63 Antibody (500X)	20μL	1tube
6	Substrate Solution	12mL	1vial
7	Stop Solution (2N H <sub>2</sub> SO <sub>4</sub> )	6mL	1vial
8	Plate Seals		2sheets

\* Sufficient to create 4 standard curves with n=2.

#### Required Materials Not Included in the Kit

- Micropipettes (10 ~ 1000  $\mu$ L)
- Multichannel micropipette
- Multichannel micropipette Reservoir
- Plate shaker
- Microplate reader (enable to measure at wavelength 450nm)
- Plate washer

---

## 【II】 Preparation of Reagents and Samples

---

### 【II – 1】 Preparation of Washing Buffer

- Dilute Washing Buffer (10 $\times$ ) to x10 with purified water.  
e.g. For 1 plate, add 225 mL of purified water to 25mL of Washing Buffer (10 x) and mix well.

### 【II – 2】 Preparation of CD63 Standard Beads solution

	Concentration (ng/mL)	CD63 Standard Beads	Assay Buffer	Dilution factor
A	400			
B	40	50 $\mu$ L of A	450 $\mu$ L	10
C	20	250 $\mu$ L of B	250 $\mu$ L	2
D	10	250 $\mu$ L of C	250 $\mu$ L	2
E	5	250 $\mu$ L of D	250 $\mu$ L	2
F	2.5	250 $\mu$ L of E	250 $\mu$ L	2
G	1.25	250 $\mu$ L of F	250 $\mu$ L	2
H	0.625	250 $\mu$ L of G	250 $\mu$ L	2

- To prepare Solution B, add 450 $\mu$ L of Assay Buffer into 50 $\mu$ L of CD63 Standard Beads (Solution A: Vortex Standard Beads for 30 seconds before use.), and then mix well (10 times dilution). To prepare Solution C, add 250 $\mu$ L of Assay Buffer into 250 $\mu$ L of Solution B,

and then mix well (2 times dilution). Similarly, 2 times dilution series for Solution D through H should be prepared.

- Use 100 $\mu$ L for measurement, using 2 wells for each solution (n=2).
- Diluted CD63 Standard Beads Solution(0.625~40 ng/mL) should be freshly prepared at each time before use.

### **【 II – 3】 Preparation of antibody solution**

- Dilute HRP conjugated anti-CD63 antibody (500x) to 500 folds using Assay Buffer.  
e.g. For 1 plate, add 20 $\mu$ L of antibody (500x) into 10mL of Assay Buffer. Mix by inverting the tube.
- \* Diluted antibody solution should be freshly prepared at each time before use.

### **【 II – 4】 Preparation of Samples (For cell culture medium supernatant)**

The culture medium containing 10% fetal bovine serum (FBS) can be directly used to quantitate exosomes. Grow cells until it reach to the confluent state. Collect the medium supernatant, centrifuge at 2000 xg for 10min, and then used as a sample. If necessary, dilute with Assay Buffer.

### **【 II – 5】 Preparation of Samples (For serum)**

Dilute the serum sample at least four times with Assay Buffer.

- \* Out of 100 $\mu$ L reaction solution, serum should not be exceeded more than 25 $\mu$ L.

### **【 II – 6】 Sample Storage**

Samples should be stored at 2~8 $^{\circ}$ C after preparation until measurement.

---

## **【 III】 Sample measurement procedure**

---

1. Bring anti-CD63 antibody solid phased plate and the reagents to the room temperature.
2. Prepare CD63 Standard Beads solution by serial dilution. (from step 【 II – 2】 )

3. Add 100 $\mu$ L each of serial diluted CD63 Standard Beads solution (0.625 ~ 40ng/mL) or Sample solution into the well.
4. Seal the microplate with Plate Seals, place into plate shaker, and then shake it at 800 rpm for 30 sec.
5. Incubate at room temperature for 2hrs for static reaction.
6. Discard all the reaction solution, and then rinse each well with 300 $\mu$ L of Washing Buffer (from step 【II-1】 ). Repeat this step for 3 times.
7. Add 100 $\mu$ L each of diluted HRP conjugated anti-CD63 antibody (from step 【II-3】 ) to the well.
8. Seal the plate, and then shake in the plate shaker.
9. Incubate at room temperature for 2hrs for static reaction.
10. Discard the antibody solution, and then rinse each well with 300 $\mu$ L of Washing Buffer. Repeat this step for 3 times.
11. Add 100 $\mu$ L of Substrate Solution into each well, and then incubate at room temperature protected from light for 20min for static reaction.
12. Visually confirm the coloring, and then add 50 $\mu$ L each of Stop Solution.
13. Place into the plate reader, and read the absorbance of each microwell on a spectrophotometer at the wavelength of 450nm.
14. Create a standard curve for the CD63 protein by plotting the mean absorbance (y axis) against the CD63 Standard concentration (x axis) (Fig. 1).


Fig.1 Example of Standard curve for CD63 Standard Beads.

## 【IV】 Example of Results

### 【IV- 1】 Cell culture medium supernatant

The culture medium of colon cancer cell line, HCT116, was collected, and exosomes were purified by ultracentrifugation. The purified exosomes, 15.6, 31.3, 62.5, 125, 250, and 500 ng/mL each, were added to the well, and measured using this kit (Fig.2).


Fig.2 Exosomes from the culture supernatant of HCT116

### 【IV- 2】 Standardization of measured values

The Standard curve was created using the measurement of CD63 Standard Beads (Fig. 3a). From here (Fig.3a), assuming 10ng/mL of CD63 Standard Beads as 1 unit/mL, the OD450 for 1unit will be 0.6. For the exosomes purified from HCT116, the protein concentration for 1 unit of OD450, which will correspond to 0.6, is about 390ng/mL (Fig.3b). Therefore, 390 ng/mL of HCT116 derived exosomes can be considered as 1unit/mL of CD63 positive exosomes.

Presenting the exosome measurement by the units, we could standardize it and/or normalize the measurements, and then will be able to compare exosome measurements directly between different samples or different experiments.


Fig.3 Standardization of measured values

---

### **【Reference】**

---

- 9) J. Skog, T. Wurdinger, S. van Rijn, D. H. Meijer, L. Gainche, M. Sena-Esteves, W. T. Jr. Curry, B. S. Carter, A. M. Krichevsky and X. O. Breakefield: *Nat Cell Biol.*, **10**, 1470 (2008).
- 10) T. Pisitkun, R. F. Shen and M. A. Knepper: *Proc Natl Acad Sci USA.*, **101**, 13368 (2004).
- 11) S. Runz, S. Keller, C. Rupp, A. Stoeck, Y. Issa, D. Koensgen, A. Mustea, J. Sehouli, G. Kristiansen and P. Altevogt: *Gynecol Oncol.*, **107**, 563 (2007).
- 12) S. Keller, A. K. Konig, F. Marme, S. Runz, S. Wolterink, D. Koensgen, A. Mustea, J. Sehouli and P. Altevogt: *Cancer Lett.*, **278**, 73 (2009).
- 13) C. Lasser, V. Seyed AlikhS. Gabrielsson, J. Lotvall and H. Valadi: *J Transl Med.*, **9**, 9 (2011).
- 14) Y. Naito, Y. Yoshioka, Y. Yamamoto and T. Ochiya: *Cell Mol Life Sci.*, **74**, 697 (2017).
- 15) A. Zoraida and M. Yáñez-Mó: *Front Immunol.*, **5**, 442 (2014).
- 16) V. Filipe, A. Hawe and W. Jiskoot: *Pharm Res.*, **27**, 796 (2010).

**COSMO BIO CO., LTD.**

Inspiration for Life Science

Licensed by Cosmo Bio Co., LTD.